



InVitro Biotainers, PET, sterile, Thermo Scientific Nalgene



55

Ideal for storage of pH sensitive samples.

- Excellent gas barrier properties
- Provides safe storage for samples from -40°C to 70°C
- Ideal for storage of media culture
- Sterile
- Moulded in polyethylene terephthalate copolyester
- Non-cytotoxic and non-pyrogenic
- FDA compliant (21CFR)
- Meet USP Class VI and European Pharmacopeia
- Modified Abnormal Toxicity Test

Catalogue No	Alt. No	Description	Material	Capacity	Inner pack qty	Pack qty
BTR-985-010M	3025-42	InVitro Biotainer	PETG	125mL	5	100
BTR-985-030G	3005-42	InVitro Biotainer, with handgrip	PETG	500mL	5	70
BTR-985-040D	3005-70	InVitro Biotainer, with handgrip	PETG	500mL	35	70
BTR-985-060U	3110-42	InVitro Biotainer, with handgrip	PETG	1L	5	35
BTR-985-050A	3110-35	InVitro Biotainer, with handgrip	PETG	1L	35	35
BTR-985-090L	3230-42	InVitro Biotainer, with handgrip	PETG	2L	5	20
BTR-985-080X	3230-20	InVitro Biotainer, with handgrip	PETG	2L	20	20
BTR-986-060X	3415-42	InVitro Biotainer, without handle	PETG	5L	1	6
BTR-986-050R	3415-16	InVitro Biotainer, with handle	PETG, polyethylene handle	5L	1	6



InVitro Biotainers, polycarbonate, sterile, Thermo Scientific Nalgene



55

- Sterile
- Strong, impact resistant
- Silicone lined polypropylene closures
- Printed graduations (except 5mL vial)
- Space saving square shape, with ribbed hand grips (except 125mL)
- Meets USP 87, 88; non-pyrogenic and non-cytotoxic
- Meets USP Class VI and European Pharmacopoeia Modified Abnormal Toxicity Test standards

Sterile and ready to use containers, moulded in blue tinted polycarbonate, that provide safe storage from -100°C to +99°C.

Catalogue No	Alt. No	Description	Material	Capacity	Inner pack qty	Pack qty
BTR-986-040U	3500-05	InVitro Biotainer	Polycarbonate	5mL	100	500
BTR-985-020J	3030-42	InVitro Biotainer with handgrip	Polycarbonate	125mL	5	50
BTR-985-070R	3120-42	InVitro Biotainer with handgrip	Polycarbonate	1L	5	35
BTR-985-100L	3233-42	InVitro Biotainer with handgrip	Polycarbonate	2L	5	20
BTR-985-120F	3405-16	InVitro Biotainer carboy with handle	Polycarbonate, polyethylene handle	5L	1	6
BTR-985-140W	3405-66	InVitro Biotainer carboy with handle	Polycarbonate	5L	6	6
BTR-985-130C	3405-42	InVitro Biotainer carboy without handle	Polycarbonate	5L	1	6
BTR-985-110Y	3405-06	InVitro Biotainer carboy without handle	Polycarbonate	5L	6	6
BTR-985-160Q	3410-42	InVitro Biotainer carboy without handle	Polycarbonate	10L	1	2
BTR-985-150T	3410-08	InVitro Biotainer carboy with handle	Polycarbonate, polyethylene handle	10L	1	2
BTR-985-170N	3423-42	InVitro Biotainer carboy without handle	Polycarbonate	10L	1	3