

Pipettes, serological, polystyrene

Fisherbrand®

S 41

Suitable for tissue culture work.

A black magnifier stripe runs lengthwise along the pipette making the meniscus obvious and reducing reading errors. The pipettes have excellent optical clarity with high contrast black markings making it easier to read volumes. The 5mL, 10mL, 25mL and 50mL sizes have both ascending and descending graduations. Tip design on the 25mL and 50mL sizes increases the effect of surface tension to minimise dripping. The individually wrapped pipettes come in a clear, all plastic wrap and the multipacks come in a recyclable LDPE bag.

Individually wrapped

Catalogue No	Capacity, mL	Graduations, mL	Colour	Inner pack qty	Pack qty
FB55481	1	0.01	Yellow	100	1,000
FB55482	2	0.01	Green	100	500
FB55483	5	0.1	Blue	50	200
FB55484	10	0.1	Orange	50	200
FB55485	25	0.2	Red	50	200
FB55486	50	0.5	Purple	25	100

Multipack

Catalogue No	Capacity, mL	Graduations, mL	Colour	Inner pack qty	Pack qty
FB55488	1	0.01	Yellow	25	1,000
FB55489	2	0.01	Green	25	500
FB55490	5	0.1	Blue	50	500
FB55491	10	0.1	Orange	25	500
FB55492	25	0.2	Red	20	200
FB55493	50	0.5	Purple	20	100



Pipettes, serological, polystyrene, plastic wrap, Stripette®

CORNING

S 66

Ideal for sterile applications such as tissue culture.

- Feature a unique vertical strip which creates a thermometer-like effect at the meniscus, making volume easy to read
- 'Banana-peel' plastic wrap provides easy, aseptic opening and pipette removal
- Individual, clear, plastic wrap allows full view of pipette for easy size identification and reduces storage space requirements
- **TKV-671-090B** has a wide tip for handling viscous fluids
- Sterile, non-pyrogenic and DNase/RNase free
- Exclusive anti-drip tip assures accurate delivery
- Bidirectional graduations provide choice of ascending and descending scales

Catalogue No	Alt. No	Volume, mL	Colour	Inner pack qty	Pack qty
TKV-671-010C	4011	1	Yellow	100	1,000
TKV-671-020W	4012	1	Yellow	100	200
TKV-671-030T	4021	2	Green	100	1,000
TKV-671-040Q	4051	5	Blue	50	200
TKV-671-050N	4101	10	Orange	50	200
TKV-671-060K	4492	10	Orange	50	200
TKV-671-070H	4251	25	Red	50	200
TKV-671-080E	4501	50	Purple	25	100
TKV-671-090B	4484	100	Aqua	10	100



Pipettes, serological, polystyrene, plastic wrap, bulk pack, Stripette®

CORNING

S 66

- Sterile, non-pyrogenic and DNase/RNase free
- Bidirectional graduations provide choice of ascending and descending scales
- Exclusive anti-drip tip assures accurate delivery

Catalogue No	Alt. No	Volume, mL	Colour	Inner pack qty	Pack qty
TKV-673-011L	4010	1	Yellow	50	1,000
TKV-673-031F	4020	2	Green	50	1,000
TKV-673-051W	4050	5	Blue	50	500
TKV-673-071Q	4100	10	Orange	50	500
TKV-673-081N	4250	25	Red	50	200
TKV-673-091K	4500	50	Purple	50	100

