



Pipettes, serological, polystyrene, paper/plastic wrap, Stripette®

CORNING



66

Ideal for sterile tissue culture applications, especially while wearing latex gloves as the wrap reduces static cling.

- Individually wrapped
- A strip creates a thermometer-like effect at the meniscus, making the volume easy to read
- Colour coding on wrapping allows easy size identification
- Sterile, non-pyrogenic and DNase and RNase free
- Exclusive anti-drip tip assures accurate delivery
- Bidirectional graduations provide choice of ascending and descending scales

Catalogue No	Alt. No	Volume, mL	Colour	Inner pack qty	Pack qty
TKV-670-011G	4485	1	Yellow	50	1,000
TKV-670-031A	4486	2	Green	50	1,000
TKV-670-051R	4487	5	Blue	50	200
TKV-670-071L	4488	10	Orange	50	200
TKV-670-091F	4489	25	Red	25	200
TKV-670-111C	4490	50	Purple	25	100
TKV-670-121W	4491	100	Aqua	10	100



Pipettes, serological, cleanroom pack, Stripette®

CORNING



66

Ideal for sterile tissue culture applications, especially while wearing latex gloves, as the paper/plastic wrap reduces static cling.

- Triple wrapped: inner 'dust cover' liner for the carton, a second sealed sleeve over each inner package, pipettes are individually wrapped
- Sterile, certified non-pyrogenic and DNase and RNase-free
- Unique colour-coded magnifying stripe for easy meniscus viewing
- Calibrated to deliver (TD-EX 20°C) contents with blow out, accuracy within $\pm 2\%$ at full volume
- Anti-drip tip reduces accidental drips and ensures accurate delivery
- Clean room pack style eliminates the need for alcohol wipe downs and similar procedures
- Colour coded pipette wraps include lot numbers and allow full view and easy opening by either pop through or peel apart

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Colour	Length, mm	Pack qty
TKV-680-010B	7041	1	0.01 (interval) 0.2 (negative)	Yellow	278	1,000
TKV-680-020V	7042	2	0.01 (interval) 0.2 (negative)	Green	277	1,000
TKV-680-030S	7045	5	0.01 (interval) 2.5 (negative)	Blue	348	200
TKV-680-040P	7015	10	0.1 (interval) 3 (negative)	Orange	344	200
TKV-680-050M	7016	25	0.2 (interval) 10 (negative)	Red	344	200
TKV-680-060J	7017	50	0.5 (interval) 10 (negative)	Purple	370	100
TKV-680-070G	7000	100	1 (interval)	Aqua	348	100

