

Cell/Tissue Culture and Cell Transfection

Apparatus and consumables - TC-treated slides


Chamber slides, TC-treated, Thermo Scientific Nunc®, Lab-Tek™ II, Chamber Slide™



82

- Treated for excellent attachment and growth of cells
- Inert hydrophobic well border printed on slide
- Removable, polystyrene chamber
- Slide separator included in each pack
- Non-fluorescent microscope slide, glass (25mm x 75mm x 1.2mm)
- Slide material: glass

Catalogue No	Alt. No	No. of wells	Volume, working, mL	Culture area, cm²	Inner pack qty	Pack qty
TKT-210-130A	154453	1	2.0 to 4.5	8.6/well	16	96
CBR-160-001Y	154453SP	1	2.0 to 4.5	8.6/well	-	16
TKT-210-150R	154461	2	1.0 to 2.0	4.0/well	16	96
CBR-160-002G	154461SP	2	1.0 to 2.0	4.0/well	16	16
TKT-210-170L	154526	4	0.5 to 1.0	1.7/well	16	96
CBR-210-170M	154526SP	4	0.5 to 1.0	1.7/well	-	16
TKT-210-190F	154534	8	0.2 to 0.5	0.7/well	16	96
CBR-160-008R	154534SP	8	0.2 to 0.5	0.7/well	-	16



Chamber slides, TC-treated, Thermo Scientific Nunc®, Lab-Tek™ II-CC2™, Chamber Slide™



82

- Chemically modified growth surface on glass slide that mimics polylysine
- Growth surface remains stable without refrigeration
- Provides binding sites optimal for fastidious growth cells (e.g. neurons)
- Light blue frosted writing area
- Slide separator included in each pack

Catalogue No	Alt. No	No. of wells	Tray quantity	Volume, working, mL	Culture area, cm²	Inner pack qty	Pack qty
TKT-215-010D	154739	1	8	2.0 to 4.5	8.6/well	16	96
CBR-161-001C	154739SP	1	8	2.0 to 4.5	8.6/well	-	16
TKT-215-020A	154852	2	8	1.0 to 2.0	4.0/well	16	96
CBR-161-002A	154852SP	2	8	1.0 to 2.0	4.0/well	-	16
TKT-215-030U	154917	4	8	0.5 to 1.0	1.7/well	16	96
CBR-161-004T	154917SP	4	8	0.5 to 1.0	1.7/well	-	16
TKT-215-040R	154941	8	8	0.2 to 0.5	0.7/well	16	96
CBR-161-008L	154941SP	8	8	0.2 to 0.5	0.7/well	-	16



Chambered coverglass, Thermo Scientific Nunc®, Lab-Tek™ II



76

For adherent cell culture.

- Excellent for confocal image analysis
- Optimal for high power inverted microscope viewing
- Borosilicate glass No. 1.5 thickness (0.16mm to 0.19mm)
- Non-removable cell culture grade polystyrene medium chamber
- Cover with clear polystyrene handle for aseptic handling

Catalogue No	Alt. No	No. of wells	Culture area, cm²	Volume, working, mL	Pack qty
CBR-177-382N	155360SP	1	8.6	2.0 to 4.5	16
TKT-230-010T	155360	1	8.6	2.0 to 4.5	96
CBR-177-406C	155379SP	2	4	1.0 to 2.0	16
TKT-230-020Q	155379	2	4	1.0 to 2.0	96
CBR-177-330J	155382SP	4	1.7	0.5 to 1.0	16
TKT-230-030N	155382	4	1.7	0.5 to 1.0	96
CBR-177-370U	155409SP	8	0.7	0.2 to 0.5	16
TKT-230-040K	155409	8	0.7	0.2 to 0.5	96



All products are for research use only, unless otherwise indicated.