

Tubes, PCR*, 0.2mL and 0.5mL, flat and domed cap



- Polypropylene single tubes with attached caps
- Compatible with leading thermal cyclers with heated lids
- Unique cap design closes and opens easily without tools yet ensures a tight fit to reduce sample evaporation
- The flat caps, in combination with the various colours, allow simple labelling and fast sample classification
- In packs of 1,000

Catalogue No	Alt. No	Volume, mL	Colour	Cap style
TUP-113-010W	781300	0.2	Clear	Domed
TUP-114-010Q	781305	0.2	Clear	Flat
TUP-114-020N	781301	0.2	Pink	Flat
TUP-114-030K	781302	0.2	Yellow	Flat
TUP-114-040H	781303	0.2	Green	Flat
TUP-114-050E	781304	0.2	Blue	Flat
TUP-115-010K	781310	0.5	Clear	Flat
TUP-115-020H	781311	0.5	Pink	Flat
TUP-115-030E	781312	0.5	Yellow	Flat
TUP-115-040B	781313	0.5	Green	Flat
TUP-115-050V	781314	0.5	Blue	Flat

LB



Tube strips, PCR*, 8 tube with cap strips, 0.2mL, domed cap



- Polypropylene 8 tube strips with attached single caps
- Compatible with leading thermal cyclers with heated lids
- Domed caps facilitate one-handed operation
- Positive seal, easy to cap and open

Catalogue No	Alt. No	Description	Colour	Volume, mL	Pack qty
TUP-117-010V	781332	Standard	Clear	0.2	120
TUP-118-010P	781333	Low profile	Clear	0.2	120

LB



Tube strip, PCR*, 8 tube, 0.2mL, domed cap



- Polypropylene 8 tube strips with attached cap strips
- Compatible with leading thermal cyclers with heated lids
- Domed caps facilitate one-handed operation
- Positive seal, easy to cap and open

Catalogue No	Alt. No	Colour	Pack qty
TUP-116-010E	781330	Clear	125

LB



*Polymerase Chain Reaction (PCR) is a process covered by patents owned by Hoffman-La Roche