

Aspirating systems, VACUBOY

INTEGRA

a 04

Enables liquid residues to be disposed of reliably, precisely and conveniently:

- With a precisely dosed volume
- Without loss of drops
- Without risk of contamination
- Without causing tiredness
- Energy saving, the manual operator - the heart of the system - can be operated as a single or 8 channel aspirator
- All parts in contact with fluids are autoclavable

VACUBOY hand operator set comprises of: VACUBOY hand operator, stand, rubber adapter for Pasteur pipettes, 1 channel stainless steel aspiration adapter 40mm, 8 channel stainless steel aspiration adapter 40mm, 1 channel plastic aspiration adapter for disposable tips, 1 channel plastic aspiration adapter for disposable tips with ejector, 8 channel plastic aspiration adapter for disposable tips with ejector. Can be used with VACUSAFE comfort vacuum system. Full range of accessories available.



Catalogue No	Alt. No	Description
PMR-150-010N	155 510	VACUBOY hand operator, comprises of hand operator, rubber adapter and 1 channel stainless steel aspiration adapter 40mm
PMR-150-035U	155 500	VACUBOY hand operator set

Accessories

Catalogue No	Alt. No	Description
PMR-150-500T	155 522	Single channel aspiration adapter 150mm, flame sterilisable
PMR-150-503N	155 520	8 channel plastic aspiration adapter for disposable tips, with ejection device
PMR-150-507F	155 503	8 channel stainless steel adapter, flame sterilisable
PMR-150-509B	155 502	Single channel stainless steel adapter 40mm, flame sterilisable
PMR-150-511X	155 504	Single channel plastic adapter for pipettor tips
PMR-150-513K	155 524	4 channel stainless steel aspiration adapter
PMR-150-515G	155 525	Single channel stainless steel aspiration adapter, flame sterilisable, 280mm
PMR-150-517C	155 505	Rubber adapter for Pasteur pipettes
PMR-150-525D	155 526	Single channel tip adapter
PMR-150-520N	155 501	Stand VACUBOY for hand operator
PMR-100-610N	155 065	PIPETBOY and VACUBOY stand

