

Microplates/strips polystyrene, 96 well, Thermo Scientific Immulon®

Thermo
SCIENTIFIC

04

Designed for use in immunoassay applications.

Immulon® 1B - exhibits a consistent medium binding affinity for hydrophobic proteins and complexes.

Immulon® 2HB - high binding, irradiated to provide increased binding affinity for hydrophilic proteins and complexes.

Immulon® 4HBX - extra high binding, irradiated, offers maximum protein uptake and guaranteed coefficient of variation (CV).

Catalogue No	Alt. No	Binding property	Well shape	Pack qty
DIS-950-010M	3355	Immulon® 1B	Flat	50
DIS-950-030G	3555	Immulon® 1B	U	50
DIS-950-050A	3455	Immulon® 2HB	Flat	50
DIS-950-070R	3655	Immulon® 2HB	U	50
DIS-950-090L	3855	Immulon® 4HBX	Flat	50
DIS-950-100L	8404	UV transparent	Flat	40

Immulon® 1B Flat

Catalogue No	Alt. No	Description	Format	Well shape	Pack qty
DIS-951-005W	6310	Removawell assembly	12 wells/strip	Flat	100
DIS-951-010G	6301	Removawell strip	12 wells/strip	Flat	320
DIS-955-010F	6505	Dividastrip assembly	2 strips/8 wells	Flat	100

Immulon® 1BC


Catalogue No	Alt. No	Description	Format	Well shape	Pack qty
DIS-953-010R	6381	Removawell assembly	8 wells/strip	C	100

Immulon® 2HB

Catalogue No	Alt. No	Description	Format	Well shape	Pack qty
DIS-951-020D	6309	Removawell assembly	12 wells/strip	Flat	100
DIS-951-030A	6302	Removawell strip	12 wells/strip	Flat	320
DIS-953-030L	6382	Removawell assembly	8 wells/strip	C	100
DIS-955-030W	6506	Dividastrip assembly	2 strips/8 wells	Flat	100

Immulon® 4HBX

Catalogue No	Alt. No	Description	Format	Well shape	Pack qty
DIS-951-080Y	6405	Removawell assembly	12 wells/strip	Flat	100
DIS-951-090F	6404	Removawell strip	12 wells/strip	Flat	320
DIS-951-100F	6484	Removawell assembly	8 wells/strip	C	100



Microplates/strips, Thermo Scientific Microlite®

Thermo
SCIENTIFIC

S

04

- Microlite® 1+ Densely pigmented white polystyrene offering medium binding with enhanced dynamic range and very low cross talk
- Microlite® 2 An irradiated plate for luminescence assays requiring higher binding, high reflectivity and minimal cross talk
- Microlite® 2+ Irradiated, offers high binding with enhanced dynamic range and very low cross talk
- Microlite® TCT A 96 well sterile plate for luminescent tissue culture assays

Enhances the sensitivity of luminescent reactions by increasing signal reflectivity. Flat well shape.

96 well microplates

Catalogue No	Alt. No	Binding property	Well shape	Pack qty
DIS-956-010W	7571	Microlite® 1+	Flat	50
DIS-956-110S	7572	Microlite® 2+	Flat	50

Removawell assemblies, 12 wells/strip

Catalogue No	Alt. No	Binding property	Well shape	Pack qty
DIS-956-030Q	7561	Microlite® 1+	Flat	100
DIS-956-130M	7562	Microlite® 2+	Flat	100

Removawell strips, 12 wells/strip

Catalogue No	Alt. No	Binding property	Well shape	Pack qty
DIS-956-040N	7566	Microlite® 1+	Flat	320
DIS-956-140J	7567	Microlite® 2+	Flat	320